

CEQUENCE

Cervical Made Simple




Product Range & Sizing

Plate Length	1 level: 10mm- 26mm (even sizes only), 2-5 levels:14mm- 22mm (length is measured from screw hole to screw hole)
Screws	Self-Drilling/Self-Tapping 4mm x 12mm-18mm , Flat Tip/Self-Tapping 4mm x 12mm-18mm, Rescue 4.5mm x 12mm-18mm (even sizes only)
Fixation Style	Semi-Constrained (20° Conical ROM) Constrained (0° Conical ROM)

Visible Locking Mechanism

Nominal Endhole Screw Angle
10° Cephad/Caudal & 10° Medial



One Driver for Locking Mechanism and Bone Screws



CEQUENCE- Features & Benefits

- Semi-constrained, constrained, and hybrid plate configurations
- 1-5 level plates in a variety of lengths
- Self-drill/self-tapping screws, flat tip/self-tapping screws, and rescue screws
- Screw diameter 4mm or 4.5mm
- One torx driver for locking mechanism and bone screws