

CIMPLICITY™


Cervical Made Simple.



SPINE
SMITH

Product Range and Sizing

Footprints (Width x Depth)	13mm (W) x 10mm (D), 17mm (W) x 14mm (D)
Lordosis	5°
Heights	5.5mm, 6.5mm, 7.5mm 8.5mm, 9.5mm or 10.5mm



CIMPLICITY- Features & Benefits

- Multiple implant holder **options**
- **Radiolucent** and biocompatible PEEK
- Radiopaque tantalum markers provide **verification of implant** positioning
- Large **graft area** to aid with fusion potential
- “**Stop**” and “**Stopless**” rasps, trials and pushers
- “**Stop**” style design straddles caspar pins for **efficient unihibited usage**