



# IN:C2™

Now Open for Fusion



## Product Range & Sizing

Footprints (Width X Length)	Small: 11.5 mm (L) x 14 mm (W) Medium: 13 mm (L) x 16 mm (W) Large: 15 mm (L) x 18 mm (W)
Lordosis	6°
Height	5, 6, 7, 8, 9, 10 or 11 mm
Plating Options	Zero Profile or Anterior Fixation
Screws	CorticoCancellous, Dual Lead or Resuce
Fixation Style	Constrained



## IN:C2™- Features & Benefits

- Modular plating options with multiple screw placement and trajectories
- **Radiolucent** and biocompatible
- Locking **snap-on** face plate
- **CorticoCancellous** screw option
- Open window **maximizes graft placement**