


# CYNCH T-PLIF


Lumbar Made Simple.



SPINE  
SMITH

## Product Range & Sizing

Footprints (Width x Length)	9mm (W) x 23mm (L), 9mm (W) x 28mm (L), 11mm (W) x 23mm (L) or 11mm (W) x 28mm (L)
Shape	Convex (0° or 6° Lordosis)
Heights	7mm, 8mm, 9mm, 10mm, 11mm, 12mm or 13mm



## CYNCH T-PLIF- Features & Benefits

- A **transverse or bilateral** approach can be used for implant placement
- **Radiolucent** and biocompatible PEEK
- Radiopaque tantalum markers provide **verification of implant** positioning
- Single implant holder
- Suture Guidance allows for the ability to “**retrieve and position**” on the implant