



CYNCH TLIF

Lumbar Made Simple.



Product Range & Sizing

Footprints (Width x Length)	11mm (W) x 25mm (L) 11mm (W) x 30mm (L)
Lordosis	5°
Heights	8mm, 9mm, 10mm, 11mm 12mm or 13mm



CYNCH TLIF- Features & Benefits

- Multiple implant holder angle **options and sites** for implant manipulation
- **Radiolucent** and biocompatible PEEK
- Radiopaque tantalum markers provide **verification of implant** positioning
- Large **graft area** to aid with fusion potential
- Suture Guidance allows for the ability to “**retrieve and position**” the implant