

DESIGNED WITH  
BIOLOGY IN MIND



# VisuALIF

*Now Open for Fusion*




SPINESMITH

## Product Range & Sizing


Footprints [Width X Length]	Small: 25mm [L] x 32mm [W] Medium: 29mm [L] x 36mm [W] Large: 29mm [L] x 42mm [W]
Lordosis	5°, 10°, 15° or 20° upon special request
Height	8mm, 10mm, 12mm, 14mm
Plate	Zero Profile or Challenging Access
Screws	5.5mm x 25, 28 or 33mm
Fixation Style	Constrained



**Zero Profile Plate**



**Challenging Access Plate**



- Modular plating options with multiple screw placement and trajectories
- Radiolucent and biocompatible
- Radiopaque tantalum markers provide verification of implant positioning
- Locking snap-on face plate
- Open window maximizes graft placement